

# iQDry™

HEARING AID DRYER

## User Manual

Please read all safety instructions before operating the device.  
Keep this manual in a safe place for future reference.

### DESCRIPTION

- + 15-minute rapid dehumidification drying system
- + One-touch switch operation design
- + Easy-to-read LCD screen
- + Digital humidity and status reminder
- + 15-minute working countdown
- + 104° F upper temperature control
- + Built-in rechargeable battery
- + Battery capacity reminder
- + USB Type-C charging compatibility
- + iQDry can be used for 7 cycles before needing to be recharged
- + Easy to clean
- + Silicone anti-slip feet for stability
- + Compact, portable design

### FUNCTION

#### 15-minute rapid dehumidification and drying system

- Intelligent temperature control
- 104° F maximum temperature
- Ultra-quiet balanced fan technology for drying
- 15-minute cycle countdown
- LCD humidity indicator
- Automatic shut-off function


#### Temperature and Humidity Display

- Built-in temperature and humidity sensor
- LCD display
- Humidity values indicated by blue, green, and orange colors

- 800mAh
- Li+ Lithium Battery**
- Runs 7 cycles before needing to be recharged



### Dehumidification Drying Test



( Test Environment: Humidity 80%, Temperature 68° F )

### COMPONENTS

1. LCD display
2. Power button
3. USB Type-C port
4. Air inlet / outlet vent
5. Silicone feet



### SPECIFICATIONS

Input: 5 V=1 A Output: 5 V=1 A  
Rated Capacity: 800 mAh 2.96 W  
Weight: 3.8 oz  
Product Dimensions: 2.8 x 3.3 x 1.8 in

### Operating Instructions

#### 1. Charging the dryer

Using the included cable, connect the USB Type-C end to the dryer access port (reference 3 on photo above).


Connect the USB Type-A end to a power source. Alternatively, use a USB-A plug adapter (not included), and plug the adapter into an electrical outlet to charge the dryer.

The LCD screen displays battery capacity and charging status.

#### 2. Using the dryer

Lift the top cover and place the hearing aids or headphones in the inner tray to be dried. Close the cover, and lightly touch the power button on the LCD screen for 2-3 seconds. (reference 2 on photos)

When the dryer is on, humidity status and other details (battery capacity, checkmark, temperature, and 15-minute countdown display) will be visible on the LCD screen. (reference 1 on photos)



#### ICON INFORMATION

<b>60-100%</b>	<b>40-59%</b>	<b>39-0%</b>
Indicates humidity inside is <b>MOIST</b>	Indicates humidity inside is <b>COMFORTABLE</b>	Indicates humidity inside is <b>DRY</b>

#### 3. End of drying cycle

The product is set to operate in cycles of 15 minutes each and will automatically shut-off when the time is up. Alternatively, touch and hold the power button for 1-2 seconds to stop operation.

### SAFETY PRECAUTIONS

Handle the iQDry dryer with care. Do not use a damaged dryer, as dropping, puncturing, burning, or applying excessive pressure to the battery may impair functionality or cause injury.

**VENTING:** Do not cover or block the air inlet / outlet vent while the dryer is in operation.

**CHARGING:** Use approved adapters and USB Type-C cables to charge or power the dryer safely during charging.

**PLACEMENT:** Place the dryer on a stable, dry surface to prevent falling and damage.

**POWER MANAGEMENT:** When unplugging, grasp the plug itself rather than the cord to safely remove it.

**CLEANING AND MAINTENANCE:** Clean the dryer with a soft cloth. Do not use abrasive cleaners or spray liquids directly onto the device to avoid surface damage. When not in use, cover the dryer to prevent dust and dirt buildup.

**WATER RESISTANCE:** Do not immerse the dryer, power cord, or plug in water or any other liquids, as this may result in electric shock and equipment damage.

**CHILD SAFETY:** Keep the dryer out of reach of children to prevent misuse and potential injury.

**PROFESSIONAL SERVICE:** Do not attempt any repairs or modifications beyond routine cleaning and maintenance. Contact the manufacturer or authorized service professional for any necessary repairs.

**DAMAGE AND REPLACEMENT:** If signs of damage appear on the power cord or plug, immediately discontinue use and obtain replacement parts (preferably the manufacturer's original model) to avoid safety hazards.

**ADDITIONAL CONSIDERATIONS:** Do not place the dryer in the direct sunlight or high-temperature environments. Do not use the dryer to dry flammable or combustible materials.

By adhering to these precautions, you can ensure the safe and effective operation of your iQDry dryer lithium-ion battery for a long time.





EAR TECHNOLOGY CORPORATION®  
1-800-327-8547 | info@eartech.com | eartech.com



## WARRANTY TERMS

---

Ear Technology Corporation warrants this product (iQDry™) against any defects that are due to faulty material or workmanship for one (1) year from the date of purchase. A purchase receipt or other proof of date of the original purchase is required for service under this agreement, otherwise the date of manufacture will be used. Ear Technology Corporation will, at its sole discretion, repair or provide a new or refurbished unit. The exchange unit is warranted for the remainder of the original warranty.

This warranty only covers failures due to defects in materials or workmanship which occur during normal use. It does not cover failures which result from accident, misuse, abuse, neglect, mishandling, misapplication, improper cleaning, line power

surge, or damage that is attributable to acts of God. Warranty is void on any unit that has been modified or tampered with.

Limitation: Ear Technology Corporation shall not be liable for incidental or consequential damages resulting from use of this product, or arising out of any breach of this warranty. All express and implied warranties, including the warranties of merchantability and fitness for a particular purpose, are limited to the applicable warranty period set forth above.

This warranty gives you specific rights, and you may also have other rights which vary by state or other government jurisdiction.

## HOW TO OBTAIN WARRANTY SERVICE

---

For more information and details regarding troubleshooting matters in North America, call Customer Service at **1-888-327-1299** or email **service@eartech.com**

If troubleshooting points to a problem that requires repair, Ear Technology Corporation will provide you with shipping instructions and a Return Merchandise Authorization number (RMA).