
# Middle Ear Analyzer

## QUICK USER'S HANDBOOK

# ME716m

CE  
0051

resonance 

## Release history

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this page is intentionally blank

## Table of contents

|   |     |
|---|-----|
| 1. Introduction and intended use .....    | 77  |
| 2. Description .....                      | 78  |
| 2.1 Precaution .....                      | 78  |
| 2.2 Warranty and maintenance .....        | 79  |
| 2.3 Unpacking and checking the R15C ..... | 80  |
| 2.4 Accessories .....                     | 80  |
| 2.5 Connecting the accessories .....      | 84  |
| 2.6 Power Sources .....                   | 85  |
| 2.7 Installation .....                    | 85  |
| 2.8 Battery indicator .....               | 85  |
| 2.9 Symbols present on device .....       | 86  |
| 3. Components and connectors .....        | 87  |
| 3.1 Power on and off the device .....     | 88  |
| 3.2 PIN Code request .....                | 88  |
| 4. Performing tests .....                 | 89  |
| 4.1 Tympanometry .....                    | 89  |
| 4.2 Stapedial Reflex test .....           | 92  |
| 4.3 Quick .....                           | 95  |
| 4.4 ETF - Intact tympanic membrane .....  | 96  |
| 4.5 Decay .....                           | 99  |
| 4.6 ETF - Patulous Eustachian Tube .....  | 100 |
| 4.7 TEOAE Screening .....                 | 101 |
| 4.8 TEOAE Diagnostic .....                | 103 |
| 4.9 DPOAE .....                           | 105 |
| 4.10 Otoscopy .....                       | 108 |
| 5. Setup settings .....                   | 111 |
| 5.1 General Setup .....                   | 111 |
| 5.2 Impedance Setup .....                 | 117 |
| 5.3 Compliance Calibration .....          | 118 |
| 5.4 OAE Setup .....                       | 121 |
| 5.5 Impedance and OAE Calibration .....   | 122 |
| 6. Patient database management .....      | 123 |
| 7. Print management .....                 | 124 |

|  |     |
|--|-----|
| 8. Data transmission to PC .....                 | 124 |
| 9. Main Labels .....                             | 126 |
| 10. Messages and warnings .....                  | 127 |
| 11. Care and maintenance.....                    | 132 |
| 11.1 Impedance Probe maintenance .....           | 133 |
| 11.2 Otoacoustic emission probe maintenance..... | 135 |
| 12. Safety regulations.....                      | 136 |
| 12.1 Electrical Safety .....                     | 136 |
| 12.2 Device handling and operation.....          | 136 |
| 12.3 Measurement Safety.....                     | 136 |
| 13. Technical specifications.....                | 137 |
| 14. Electromagnetic compatibility.....           | 141 |

# 1.Introduction and intended use

RESONANCE thanks you for having purchased one of its medical devices.

This manual contains information and warnings aimed at ensuring the correct and safe use of the R16M Middle Ear Analyzer available in the following versions:

| Description | Product code |
|-------------|--------------|
| R16M        | 4300104100   |
|             | 4300104145   |
| R16M-D      | 4300104110   |
| R16M OAE    | 4300104115   |
|             | 4300104125   |
|             | 4300104135   |
| R16M-D OAE  | 4300104120   |
|             | 4300104130   |
|             | 4300104140   |

Carefully store this manual away in a safe place for each future consultation.

The R16M Middle Ear Analyzer is designed and manufactured to meet or exceed all quality and safety requirements, and has been certified with CE-symbol according to European Medical De-vices Directive. All the functions are software-controlled, upgrading the software at a later date will be simple and free of charge.

The R16M Middle Ear Analyzer is intended to be used by audiologists, hearing healthcare professionals, or trained technicians.

A middle ear analyzer is a medical device used to assess the condition and function of the outer and middle ear. This is performed by inserting a small probe into the outer ear canal whilst presenting and measuring changes in level via a small microphone. Changes in tone level detected in the ear canal are monitored whilst the air pressure level is varied. These changes in level can be used to make inferences about the status and condition of the tympanic membrane and middle ear ossicular chain. Additional acoustic stimuli may also be presented in order to elicit reflex responses that can be of further diagnostic value.

The Otoacoustic Emissions (OAE) Screener is a device testing the cochlear functionality which is connected with the amplification activity by the external hair cells. When the cochlea is stimulated by a sound the external hair cells vibrate producing in their turn an almost imperceptible sound which goes back to the middle-ear. This subtle sound is measured by a tiny probe inserted into the ear canal. The user must therefore have an adequate knowledge of audiological diagnostic equipment, expertise on the examination methods used and knowledge of English.

This device is intended to be used in medical environments with low levels of background noise. The operator must exercise caution in interpreting results where background noise has not been sufficiently isolated.

It is recommended that the device be operated within an ambient temperature range of 15 -35 degree Celsius.

This operating manual is not intended as a training manual for audiometry and tympanometry. The reader should consult standard audiology texts for the theory and application of the screening tests provided by this device.

## 2. Description

### 2.1 Precaution

**WARNING** indicates a hazardous situation which, if not avoided, could result in death or serious injury.

**NOTICE** is used to address practices not related to personal injury.

This instrument has been built and tested in compliance with product standards and left the production factory in perfect technical safety conditions.

- Prior to testing, the performance of the unit should be checked by the operator on their own ear. This is to ensure that the unit produces results expected results.
- Ensure that the correct tip size is selected to give a good seal in the ear canal. Use only the single use ear tips provided with the unit. Replacement tips are available from RESONANCE or one of their agents.
- The device transducers and parts in contact with the patient should be cleaned following every use. Appropriate antiseptic solution should be used in compliance to the maintenance and cleaning section of this manual.
- The device must be used only with the transducers provided by the manufacturer. Use of alternative transducers, or their replacement, will require re-calibration.
- The device is intended to be powered from Li-ion internal rechargeable battery or from a medical CE power supply. Use of power supplies other than those provided by the manufacturer is not permitted. When used, position so that the power plug can be easily disconnected.
- Although the instrument fulfils the relevant EMC requirements, precautions should be taken to avoid unnecessary exposure to strong electro-magnetic fields, e.g. from mobile phones, microwaves, radiotherapy devices, etc.; these may induce malfunctions in their performance. If the device is used adjacent to other equipment it must be observed that no mutual disturbance appears.
- R16M does not have specific ingress protection. Care should be taken to avoid liquid ingress. Immersion of the device may result in serious damage.

Any possible modification and/or alteration to this device is prohibited, without prior authorization thereto by the Manufacturer.

If the R16M is connected to other devices with their own power source, the system conforms to IEC 60601-1 must be checked under the responsibility of the installer; in the event that isolation devices are used to achieve such conformity, these must comply with IEC 60601-1. Any person who connects external equipment to the device connectors has formed a system and is therefore responsible for the system to comply with the requirements of IEC 60601-1.

If in doubt, contact your RESONANCE technician or local distributor for help.

Device not to be used in areas in which flammable anaesthetics or flammable agents for disinfection or skin cleaning are applied.

Electric and electronic waste may contain hazardous substances and therefore must be disposed of separately. Such products will be marked with the crossed-out wheeled bin shown. The cooperation of the user is important in order to ensure a high level of reuse and recycling of electric and electronic waste. Failure to recycle such waste products in an appropriate way may endanger the environment and consequently the health of human beings.

## 2.2 Warranty and maintenance

RESONANCE warrants that:

- The R16M is free from defects in material and workmanship under normal use and service for a period of 36 months from the date of device's installation by authorized personnel to the first purchaser.
- Accessories are free from defects in material and workmanship under normal use and service for a period of 12 months from the date of device's installation by authorized personnel to the first purchaser.

Furthermore, this warranty shall not apply to, and M.R.S. shall not be responsible for, any loss arising in connection with the purchase or use of any RESONANCE device that has been:

- repaired by anyone other than an authorized RESONANCE service representative;
- altered in any way so as, in RESONANCE judgment, to affect its stability or reliability;
- subject to misuse or negligence or accident, or which has had the serial number altered, effaced or removed; or improperly maintained; or used in any manner other than in accordance with the instructions furnished by RESONANCE.
- This warranty is in lieu of all other warranties, express or implied, and of all other obligations or liabilities of M.R.S., and M.R.S. does not give or grant, directly or indirectly, the authority to any representative or other person to assume on behalf of M.R.S. any other liability in connection with the sale of RESONANCE products.
- The R16M may be repaired only by your dealer or by a service company authorized by RESONANCE, recommended by your dealer.
- Only original RESONANCE spare or consumable parts must be used for service or repair the R16M.

In order to ensure that your device works correctly, the R16M Middle Ear Analyzer should be checked and calibrated at least once a year. We recommend to send the R16M to your RESONANCE dealer for this service, including the probe and all of the accessories. The verification of instrument calibration must be conducted by qualified service personnel authorized by the manufacturer via appropriate instrumentation. This must provide the level control and the frequency of all signals emitted by the transducers supplied with the instrument.

Should transducers require replacement, use only parts provided by the manufacturer. Replacement of transducers will require re-calibration of the instrument.

## 2.3 Unpacking and checking the R16M

On checking for packaging and content damage immediately after the device's reception; should any sign of damage or tampering be noted, please notify the carrier immediately.

If the device has been damaged during transportation, the device should be checked for any electrical or mechanical defects.

Should the device be dropped or fall, it may sustain hidden damage that may result in hazardous operation.

In case of any doubt, please do not use the device and either contact your Resonance Dealer, or e-mail the Resonance Customer Support Service directly at [support@resonance-audiology.com](mailto:support@resonance-audiology.com)

Keep all original packaging for future use.

The identification label on the rear plate holds the serial number. This should be checked and written down for later service claims.

## 2.4 Accessories

### List of accessories

Impedance probe 226 Hz/1000 Hz  
Otoacoustics emission probe (OAE version only not available for 4300104145 code)  
Insert for Contra use (R16M only)  
IP30 insert earphones for Contra use (R16M-D only)  
Calibration cavity  
Kit Tips for impedance probe, screening type

Kit Tips for OAE probe (OAE version only)  
Medical CE approved power supply  
Cradle (optional for 4300104145 code)  
Multilingual CE Operation manual  
MDS Software  
USB pen drive including: DHR documents, MDS software installer, Multilingual CE Operation manual  
Carrying bag

### Available tests

#### Tympanometry

Tympanometry 226 Hz  
Tympanometry 100 Hz  
Ipsilateral Reflexes  
Contra Reflexes (not available for 4300104145 code)  
Quick  
\*Decay  
\*Decay Free Run  
\*ETF Intact and Patulous

### Optionals available

IP30 insert earphones for Contra use (R15C only not available for 4300104145 code)

External USB Keyboard  
Wi-Fi Thermal printer  
Kit Tips MS clinical type  
Otoscope (not available for 4300104145 code)

### Otoacoustic emission

\*\*TEOAE screening  
\*\*TEOAE diagnostic  
\*\*DPOA

### \*\*Otoscopy

\*R16M-D only  
\*\* = optional test, depending on device model

## R16M and R16M-D Accessories

| Image   | Description   | Standard /Optional                 | Plug   |
|---|---|------------------------------------|--|
|    | Probe 226 Hz/1000 Hz*   | •                                  | PROBE (3)  |
|    | Insert for Contra*  | • R16M<br>– R16M-D<br>- 4300104145 | CONTRA (4)   |
|    | IP30 insert earphones for Contra use*   | • R16M-D<br>○ R16M<br>- 4300104145 | CONTRA (4)   |
|    | Calibration cavity  | •                                  | /  |
|    | Kit- Tips screening type  | •                                  | /  |
|    | Medical CE approved power supply:<br>Model: MENB1010A0503 (xyz)<br>(xyz) = B01 changeable direct plug-in<br>(xyz) = Q01 with power cord<br>(xyz) = N01 with power cord                | •                                  | 5 VDC (1)  |
|    | As alternative:<br>Medical CE approved power supply:<br>Model: ME10A0503 (xyz)<br>(xyz) = B01 changeable direct plug-in<br>(xyz) = Q01 with power cord<br>(xyz) = N01 with power cord | •                                  | 5 VDC (1)  |
|    | Charging Cradle   | ○ 4300104145                       | /  |
|  | USB pen drive   | •                                  | For PC use   |
|  | Carrying bag  | •                                  | /  |
|  | Wi-Fi Thermal printer   | ○                                  | /  |
|  | External USB Keyboard   | ○                                  | <br>(2) |

|   |                                  |                   |  |
|---|----------------------------------|-------------------|--|
|   | Impedance Kit Tips clinical type | ○                 |  |
|  | RE-VO Video Otoscope             | ○<br>- 4300104145 | <br>(2) |

● = standard      ○ = optional      - = not available      \* = applied part

The applied parts are represented by the transducers and pushbutton which are placed in contact with the patient. They are indicated by an asterisk (\*).

### R16M OAE and R16M-D OAE Accessories

| Image   | Description  | Standard /Optional | Plug           |
|---|--|--------------------|----------------|
|    | Probe 226 Hz/1000 Hz*  | ●                  | BC/PROBE (3)   |
|    | OAE Probe*   | ●                  | BC/PROBR (3)   |
|    | Insert for Contra*   | ● R15C<br>- R15C-T | ACR/CONTRA (4) |
|    | IP30 insert earphones for Contra use*  | ● R15C-T<br>○ R15C | ACR/CONTRA (4) |
|    | Calibration cavity   | ●                  | /              |
|   | Kit Tips screening type  | ●                  | /              |
|  | Kit Tips OAE   | ●                  | /              |
|  | Medical CE approved power supply:<br>Model: MENB1010A0503 (xyz)<br>(xyz) = B01 changeable direct plug-in<br>(xyz) = Q01 with power cord<br>(xyz) = N01 with power cord | ●                  | 5 VDC (1)      |

|   |   |   |  |
|---|---|---|--|
|   | As alternative:<br>Medical CE approved power supply:<br>Model: ME10A0503 (xyz)<br>(xyz) = B01 changeable direct plug-in<br>(xyz) = Q01 with power cord<br>(xyz) = N01 with power cord | • | 5 VDC (1)  |
|  | Charging Cradle   | • | /  |
|  | USB pen drive   | • | For PC use   |
|  | Carrying bag  | • | /  |
|  | Wi-Fi Thermal printer   | ○ | /  |
|  | External USB Keyboard   | ○ | <br>(2) |
|  | Impedance Kit Tips clinical type  | ○ |  |
|  | RE-VO Video Otoscope  | ○ | <br>(2) |

● = standard      ○ = optional      - = not available      \* = applied part

The applied parts are represented by the transducers and pushbutton which are placed in contact with the patient. They are indicated by an asterisk (\*).

## 25 Connecting the accessories

Please connect all the cables before powering on the device

| Reference | Symbol  | Function   |
|-----------|---|--|
| 1         | 5 VDC   | Socket for external power supply   |
| 2         |    | USB Host port for connection of pen drive (for software update) or external keyboard, Wi-Fi dongle, otoscope |
| 3         | BC/Probe  | Socket for: Impedance probe (226 Hz/1000 Hz), OAE probe, Bone vibrator. Air pressure                         |
| 4         | Contra  | Socket for Contralateral stimulation transducer connection   |
| 5         |  | Patient response push bottom mini usb plug. This connector can be used for wired data transfer               |
| 6         |   | Not in use   |

## 2.6 Power Sources

The device can operate with its own internal rechargeable battery (Li-ion 3.6 V / 3 Ah), or connected to the medical power supply, or placed on the charging cradle.

The average recharge time of an empty battery is less than 4 hours, both with direct feeder link or with its charging cradle. Two hours are enough to reach 80% charge.

To increase the battery duration, the display brightness is reduced to 20% after 60 seconds of inactivity (no action on the touch screen).

After 2 minutes of inactivity starting from the reduction of display brightness, the device switches to sleep mode. Simply touch the screen to reactivate the device. After 8 minutes of inactivity, and if not connected to the battery charger, the instrument switches off automatically. Press the ON / OFF button to restart the device.

## 2.7 Installation

Ensure all transducers are connected prior to starting the device.

Turn the device on and check the battery status: at first switching on, the battery indicator may provide an approximate value. A full charge cycle is suggested prior to its use.

## 2.8 Battery indicator

The battery charge status is showed on the display with a dedicated icon. When less than 10%, the battery indicator flashes to underline that recharging is necessary.

At low charge levels, the device displays a Low Battery warning message. Some tests may not be available following the low-battery warning.

Although the device is able to operate even when the battery is running very low, it is not recommended to use it with less than 5% battery autonomy and it is recommended to connect it to the medical power supply or to its charging cradle.

## 2.9 Symbols present on device

|   |   |
|---|---|
|  | 5 VDC Direct current (Medical CE approved power supply)   |
|  | Refer to instruction  |
|  | Refer to manual for transducer and power supply connection  |
|  | Class II device<br>(Symbol refers to external medical power supply with double insulation)                  |
|  | BF type applied parts   |
|  | Electric and electronic waste may contain hazardous substances and therefore has to be collected separately |
|  | Together with name and address indicates the manufacturer of the instrument                                 |
|  | Conformity marks to directive 93/42/CEE – notification body 0051 (IMQ)                                      |
|  | Year of production of the device  |
|  | USB Host port   |
| PAT. RESP.  | Patient response push button plug   |

### 3. Components and connectors

|   | Description  |
|---|--|
| 1 | Socket for Probe: probe 226 Hz/1000 Hz/OAE Probe + air pressure tube                                     |
| 2 | Socket for Contralateral stimulation transducer connection   |
| 3 | 7" TFT Resistive touch screen display  |
| 4 | External body  |
| 5 | Main switch  |
| 6 | Led: shows that battery is charging  |
| 7 | Socket for external medical power supply   |
| 8 | Not in use   |
| 9 | USB Host port for connection of pen drive (for software) or external keyboard. Wi-Fi dongle and otoscope |

## 3.1 Power on and off the device

To turn on the device, press briefly the Main switch (6) and wait for the Pin code request on the display. To switch off the device, press the Main switch (6) for about 2 "and then release the button.

Before powering on the device, ensure that all the accessories are already plugged-in.

## 3.2 PIN Code request

When switching the device on a PIN code is requested to access both the device functionalities and the data it contains. The default PIN code set by the manufacturer is 1234. After the PIN code has been typed in the starting screen is presented. The PIN code can be changed or disabled through the dedicated function available in the setup menu. For additional details see chapter 5.1

In case wrong codes are typed in, after 3 failures the possibility of new attempts is prevented and unlocking is only possible with the PUK code provided along with the device. Should the PUK code be lost, it can be retrieved by communicating to the manufacturer the device serial and product numbers to be found on the label on the lower side of the unit.

## 4. Performing tests

### 4.1 Tympanometry

Tympanometry is used to measure the acoustic admittance (which is also known as compliance) of the tympanic membrane and middle ear at a fixed frequency over a range of pressures.

Before proceeding with the examination, it is advisable to perform otoscopy in order to check the accessibility of the external ear canal and to evaluate the condition of the tympanic membrane.


Once otoscopy has been performed, insert the probe into the external auditory canal, by choosing an appropriately sized eartip, depending on the diameter of the ear canal size.

Exercise therefore appropriate rotation and pressure on the external auditory canal to allow the probe to obtain a perfect air-seal.

The R15C measures the admittance of the tympanic membrane and middle ear by playing a continuous tone into the ear canal at either 226 Hz or 1000 Hz. The level of this tone is calibrated at 85 dB SPL (226 Hz) or 75 dB SPL (1000 Hz) into a 2 ml cavity. The sound level this produces in the ear canal is measured using a microphone and the admittance calculated from the result. In line with normal audiometric practice admittance is displayed as an equivalent volume of air in ml (for 226 Hz) or mmho/m $\bar{U}$  (for 1000 Hz).

To record the tympanogram the admittance is measured while the air pressure in the ear canal is varied by means of a small pump. The admittance peaks when the air pressure is the same on both sides of the tympanic membrane. The change of admittance with pressure is displayed graphically and is called a tympanogram.

#### Test execution


Turn on the device and press the button **TYMPANOMETRY** (see adjacent image) on the main page to enter the Impedance test.

## 1. Test parameters

Directly press on the area of the tympanogram display to change the ear to be tested (R / L). The upper part of the tympanometric path is red (R) or blue (L) colour based on the selected ear to be tested.

Press the **226 Hz** icon (see adjacent image) to select the frequency at which to run the test. By default, the device is set to 226 Hz. The available choices are two: 226 or 1000 Hz (High Frequency Tympanometry).

Two test options are available: one single tympanogram, or three in a sequence. If you choose the first option, every time you start the test the unit performs one single tympanogram, whereas if you select the other one three consecutive tympanograms are performed, and you can afterwards select the one you wish to save in the database.

The generation of the pressure variation rate (expressed in daPa/sec) into the ear canal can be set on two modes. Select **AUTO** (see adjacent image) to initially generate the pressure variation rate at maximum speed available, then reduced when approaching the peak compliance to improve its definition. Select **FAST** to generate the pressure variation rate constantly at maximum speed available of 200 daPa/sec.

Press the **+200 -400** icon (see adjacent image) to choose the start and stop pressure range. Default setting is + 200 / -400 range alternatively, +100/-200 daPa option is available.

Press the **START** button (see adjacent image) to start the exam.

At the end of the test execution, you can repeat the test or change the ear to be tested to complete the session.

The table in the centre of the display gives examination result values. These are displayed in red for the right ear, and blue for the left.

| AUTO | Rate      | AUTO |
|------|-----------|------|
| 1.12 | ECV ml    | 1.15 |
| 0.38 | GRAD      | 0.26 |
| 57   | GRAD daPa | 81   |
| 1.19 | Peak ml   | 0.81 |
| 9    | Peak daPa | 9    |

Press **ZOOM +** or **ZOOM -** icon (see adjacent image) to get a zoom of the examination graphic display. By default, the ordinate shows mx 2ml. Zoom of x0.5 / x2.5 / x3.5 / Absolute are available.

To store the test results and return to the main page, simply press the **BACK**-button icon (see adjacent image). The device automatically asks if you want to save and/or print the test out on the thermal printer.

If save is selected, the R16M will open the patient database in order that demographics can be completed.

## 2. Various indicators

This section indicates the status of the probe. Probe open, Probe closed, ready or autorun will be displayed.

Indicates Real-time Compliance value read by the probe

Indicates Real-time Pressure value read by the probe

Indicates the tympanogram type identification based on Jerger's classification automatically detected by the device and displayed in the upper right of the tympanogram area performed at 226 Hz.

### 3. Children function

When performing tympanometric tests in children you can select the “Children” function. This option distracts them from the test by showing images and animations on the display.

It is possible to switch from adult to child mode and vice versa by pressing the button you can see in the picture. Animations automatically start by pressing the start button or, if the autorun function has been activated, by inserting the probe into the ear. In the latter case it is possible to manually start the video by holding down the child button for a few seconds. In order to visualize the test result the output key has to be pressed for 6 seconds.

## 4.2 Stapedial Reflex test

It is also possible to investigate whether a stapedial reflex is present. In this case, the 226 Hz tone is used to measure the admittance of the ear, while a short tone at a different frequency is presented (the reflex stimulus). The level of this stimulus is increased in steps until the stapedial muscles respond causing the tympanic membrane to become stiffer, or a preset maximum level is reached. When the change in admittance exceeds a predetermined threshold, this constitutes a reflex and the change in admittance at that level when the stimulus is applied is displayed as a plot against time.

The stapedial reflex is measured at the static ear canal pressure that produces the maximal admittance, so reflex measurements are taken after the tympanogram is measured when the peak admittance pressure has been established.

The reflex stimulus may be produced in the ear being measured (ipsilateral mode), the opposite ear (contralateral mode) or in both ears (ipsilateral mode followed by contralateral mode). For contralateral stimulation the reflex tone is produced in a separate transducer supplied with the instrument.

### Test execution


Press the **REFLEX** icon (see adjacent image) from the main page to enter the Reflex test.

## 1. Test parameters

Directly press on the area of the Reflexes visualization to change the ear to be tested (R/L). The upper part of the track is red or blue colour based on the selected ear.

Some menu may differ for the 4300104145 code, because they do not have the contralateral transducer and some reflex stimulus frequencies.

Press the **IPSI** icon (see adjacent image) to select whether to perform Reflexes test in **IPSI** (Ipsilateral) or **CONTRA** (Contralateral) mode. The frequencies battery displayed on the upper area is reserved for **IPSI** selection, while the lower is reserved for **CONTRA**. These tests are manually performed. If you wish to perform an ipsilateral + contralateral test sequence you can set I + C.

Press the **SETUP** icon (see picture aside) to enter the Reflexes test **SETUP** page.

Enable/disable the frequencies to be tested, the start and the maximum desired stimulus presentation level (dB HL) and the duration of the same.

By setting the starting and ending levels at the same value, the test will be performed at that specific value only, as if a manual reflex search was performed.

By setting two different values, the test will be performed looking for the threshold between the two selected values.

The device provides for different maximum output levels depending on the frequency and on the Ipsilateral or Contra output; if you wish to investigate the reflex up to the available maximum levels you only need to press the Up-arrow key until setting MAX. For additional details on the maximum output levels refer to section 13 Technical Data.

Fixed level at 90 dB

Threshold search between 80 and 100dB


Threshold search between 80 dB and the maximum available level

Press the **SETUP** icon (see adjacent image) to return to the test execution page.

Press the **START** icon (see adjacent image) to start the exam.

After the test finishes, it is possible to change sides, run a contralateral test, repeat or end the session. The results can be viewed on a larger scale by pressing and holding a graph for 2 seconds. If selected in the setup menu, the slope of the rising and falling edge of the reflex in ml/sec will be displayed, too.

Pressing the **EXIT** button, a dialog box will pop up and ask if you wish to save and/or print the test out on the thermal printer.

Select the appropriate button. The R16M will show the database from where you can select the patient to assign the result to.

## 2. Various indicators

**PROBE OPEN**  
status

Indicates probe status. The icon will read either: probe open; probe closed; ready; autorun.

**8.43**  
ml

Indicates Real-time Compliance value read by the probe

**17**  
daPa

Indicates the Real-time Pressure value read by the probe.

### 3. Children function

When performing tympanometric tests in children you can select the **CHILDREN** function. This option distracts them from the test by showing images and animations on the display.

It is possible to switch from adult to child mode and vice versa by pressing the button you can see in the picture. Animations automatically start by pressing the start button or, if the autorun function has been activated, by inserting the probe into the ear. In the latter case it is possible to manually start the video by holding down the child button for a few seconds. In order to visualize the test result, the output key has to be pressed for 6 seconds.

## 4.3 Quick

R16M sequential tests are designed exclusively for rapid screening.

When pressing **START**, the device performs in rapid sequence a tympanometric test and a reflex test.

### Test execution

Switch on the device and press the **QUICK** icon (see adjacent image) on main menu to enter Quick test.


Reflex screening can include all of the settings of the full test. The user can define the parameters.

Up to three sets of parameters can be stored. To save the settings press one of the number icon on the right side of the screen (for example icon 2) for some seconds until it turns green.

Please, refer to the previous sections of this manual for icons/buttons functions.

Some menu may differ for the 4300104145 code, because they do not have the contralateral transducer and some reflex stimulus frequencies.

Indicates the probe status: probe open; probe closed; ready or autorun.

## 4.4 ETF - Intact tympanic membrane

ETF stands for "Eustachian Tube Function", this test assesses its function in patients with intact tympanic membrane.

The automatic test, described by P. Williams (Ann Otol 84: 1975, 339-343), is performed using 226 Hz probe tone, this is the most significant frequency for Tube function measure.

Test execution

Switch on the device and press the **ETF** icon (see adjacent image) on main menu to enter the ETF test.



The screen shows a single chart with a colored stripe in the diagram upper section (red+R = right, blue+L= left) to indicate the examined ear.

Press on the colored stripe to change sides.

The screen also shows the results table, zoom icon, pump speed selection and the start/stop button.

### 1. Test parameters

Directly press on the area of the tympanogram display to change the ear to be tested (R / L). The upper part of the tympanometric path is red (R) or blue (L) color based on the selected ear to be tested.

Pressing the **AUTO** icon toggles between automatic mode, slowing the pump as it passes the tympanometric peak pressure; or running as fast as possible.

The pressure range icon is used to choose either +200/-400 daPa or +100/-200 daPa

Press the **START** button (see adjacent image) to start the test.

The ETF test requires patient's cooperation. The test performs 3 tympanograms at 226 Hz, between tympanograms a message will prompt the operator to instruct the patient prior to the next tympanogram.

At the end of the test the chart shows 3 paths, named 1, 2 and 3. The first tympanogram is executed at 0 pressure as in the classic tympanometric exam, the second at the maximum positive pressure (+200 daPa); the third at the maximum negative pressure (-400 daPa).

### The test result shows:

- Diagnostic and prognostic information on the efficiency of the eustachian tube.
- Early diagnosis of middle ear disease.
- Follow-up of acute and sub-acute inflammatory conditions during therapy.
- Follow-up in miringo-plastics with transtympanic drainage to decide the removal of the drainage tubes when the tube function is returned normal.
- Correct prognosis of tympanoplasty in patients with chronic middle otitis.

The technical data used during the test are shown in the following table displayed on screen.

| Rate | ECV  |      |
|------|------|------|
| AUTO | 1.06 |      |
| daPa | Peak | ml   |
| 14   | 1    | 0.97 |
| 15   | 2    | 0.97 |
| 17   | 3    | 1.01 |

Pressing the zoom icons alters the display scale.

Pressing the **EXIT** button, a dialog box will pop up and ask if you wish to save and/or print the test out on the thermal printer.

Select the appropriate button. The R16M will show the database from where you can select the patient to assign the result to.

## 4.5 Decay

The tone decay test, also known as the Anderson test, can only be performed in patients with confirmed stapedal reflexes. A constant sound is presented to the ear at a level sufficient to elicit a reflex for a period of 10 s. The initial size of the reflex is measured. If this level is reduced by 50% within the 10 s period it is indicative of a retrocochlear pathology.

### Test execution

Switch on the device, then press the **DECAY** icon (see adjacent image) on the main page to enter the Decay test.


The test can be performed in **IPSI** or **CONTRA** mode with a 10 or 20 seconds time setting. Pressing the chart allows you to change test ear.

Press on this icon to select the **IPSI** or **CONTRA** test mode

Press on this icon to select the duration, 10 or 20 seconds

Arrows are used to adjust the stimulus level

Arrows are used to adjust the stimulus frequency

Press this icon to start the test

## 4.6 ETF - Patulous Eustachian Tube

The Decay test is used to test patients with patulous eustachian tube. The test can last 10 or 20 sec, is performed at 0 daPa and with no acoustic stimulation.

In order to perform the test, the patient has to be instructed to breathe deeply through the nose and without opening the mouth.

In case the patient is actually affected by this pathology, the graph will result as oscillating in relation to the patient's breathing.

### Test execution

Switch on the device, then press the **DECAY** icon (see adjacent image) on the main page to enter the ETF test - Patulous Eustachian Tube.

The test can be performed in either **IPSI** or **CONTRA** modes. Pressing the chart allows you to select the test ear.

Because no acoustic stimulus is presented, the ipsi selection is not relevant

Press on this icon to select the duration, 10 or 20 seconds

Press the down arrow until the stimulus level is OFF

Because no acoustic stimulus is presented, the frequency is not relevant

Press this icon to start the test

## 4.7 TEOAE Screening

The Otoacoustic Emissions Test allows to evaluate the cochlear functionality which is connected with the amplification activity by the external hair cells. The test is performed using a small probe which is inserted into the external ear canal, stimulating the cochlea by a wideband signal and measuring the response signal produced by the external cells.

The measurement is based on a statistical evaluation of the detected signals, automatically providing a result which is expressed as either “pass” or “refer”.

### Test performance

Switch the device on and press **TEOAE-S** (see picture beside) on the main page to open the test function.

The main screen displays a graph and an indicator showing the exam progress.

Press **R** (right) or **L** (left) to switch the ear to be tested.

Press **START** to start the test.

Press **PAUSE** to interrupt and resume testing.

The round indicator below the graph shows the test progress.

When the test is finalized, the results are shown on the screen:

**PASS:** test passed, valid OAE response detected.

**REFER:** test failed, no OAE response detected.

**NOT COMPLETED:** test stopped by the operator.

**Stability:** this is an indication of the stimulus signal stability. The value indicates the percentage of valid stimuli on the total.

**Artefact:** this is an indication of the signal quality assessed during the test. This percentage is calculated on the number of acquisitions during which noise has been higher than the setup threshold level.

**Noise:** test average noise.

The noise, progress and sweep indicators are updated in real time while the test is in progress.

Press **SETUP** (see picture beside) to open the test setup page.

The menu allows for the setting of a few test parameters:

Adjustable parameters:

**Max Sweep**

Max test time. By touching this field, the maximum test time, between 30 and 75 seconds, can be selected. When a valid result is obtained the test is automatically stopped.

**Artefact Level**

In this field the maximum acceptable noise level can be selected. When the noise exceeds the set value, the measurement is rejected (artefact). By touching this field, a value between 20 and 60 dB SPL can be selected.

**Stimulus Level**

Allows for the setting of the stimulus level. By touching this field, a value between 70 and 90 dB peSPL can be selected.

A maximum of 3 protocols of predefined values can be saved by keeping the numbered fields on the right side of the screen pressed for a few seconds and until they turn green.


## 4.8 TEOAE Diagnostic

In addition to the “pass” or “refer” result, the diagnostic TEOAE test provides information on the intensity and the frequency where a response has been detected.

### Test performance

Switch the device on and press **TEOAE-D** (see picture beside) on the main page to open the test function.

The main screen displays three graphs: on the left side a histogram showing the signal/noise ratio for the enabled frequencies, on the upper right the two A and B buffers and on the lower right the signal and noise response.


Press **R** (right) or **L** (left) to switch the ear to be tested.

Press **START** to start the test.

Press **PAUSE** to interrupt and resume testing.

All results are displayed on the screen when the test is finalized:

**PASS:** test passed the measurement met the set criteria. Detection criteria can be customized in the setup menu, based on the specific screening or diagnostic needs.

**REFER:** test failed, no OAE response detected. The measurement did not meet the set criteria.

**NOT COMPLETED:** test stopped by the operator.

Below the histogram bars the SNR values for the different frequencies are shown.

**Stability:** this is an indication of the stimulus signal stability. The value indicates the percentage of valid stimuli on the total.


**Artefact:** this is an indication of the signal quality assessed during the test. This percentage is calculated on the number of acquisitions during which noise has been higher than the setup threshold level.

**Noise:** test average noise.

The noise, progress and sweep indicators are updated in real time while the test is in progress.

Press **SETUP** (see picture beside) to open the test setup page.

The menu allows for the setting of a few test parameters:

#### **Adjustable parameters:**

##### **Max Sweep**

Max test time. By touching this field, the maximum test time, between 30 and 75 seconds, can be selected. When a valid result is obtained the test is automatically stopped.

##### **Artefact Level**

In this field the maximum acceptable noise level can be selected. When the noise exceeds the set value, the measurement is rejected (artefact). By touching this field, a value between 20 and 60 dB SPL can be selected.

##### **Stimulus Level**

Allows for the setting of the stimulus level. By touching this field, a value between 70 and 90 dB peSPL can be selected.

##### **Number of frequencies per pass**

Minimum number of frequency bands in which the signal to noise ratio needs to exceed the predefined maximum noise level for the test to be passed.

##### **Frequencies and SNR**

The value of the signal to noise ratio for each frequency can be set between 3, 6 or 9 dB, or the frequency band can be disabled by setting "-".

### Correlation

It is possible to define the minimum acceptable reproducibility value, or the correlation between the signals of the A and B buffers. The value can be set between 50 and 90% in steps of 10 or OFF in case you do not wish to consider this parameter.

A maximum of 3 protocols of predefined values can be saved by keeping the numbered fields on the right side of the screen pressed for a few seconds and until they turn green

## 4.9 DPOAE

The DPOAE (Distortion Product Otoacoustic Emissions) test is performed by simultaneously presenting two tone stimuli to the aim of evoking non-linear cochlear responses generating a new frequency. The variations of the two primary frequencies kept at a constant ratio allow to investigate the cochlear functionality and to create a graph, similar to an audiometry, called DP-gram.


Otoacoustic emissions prove a good mechanical functionality of the outer hair cell synapses and of the related structures, however, the information they can provide with reference to the functionality of the internal hair cells is less precise.

DPOAEs can be detected in individuals with hearing threshold up to 35-10 dB HL.

### Test performance

Switch the device on and press **DPOAE** (see picture beside) on the main page to open the test function.

The main screen displays a DP-gram where the x-axis represents the reference frequency ( $f_2$ ) and the y-axis the detected OAE level expressed in dB.


Press **R** (right) or **L** (left) to switch the ear to be tested.

Press **START** to start the test.

Press **PAUSE** to interrupt and resume testing.

The test results are displayed on the screen when the test is finalized:

**PASS:** test passed; the measurement met the set criteria. Detection criteria can be customized in the setup menu, based on the specific screening or diagnostic needs.

**REFER:** test failed, no OAE response detected. The measurement did not meet the set criteria.

**NOT COMPLETED:** test stopped by the operator.


Otoacoustic emissions are represented by the green dots and noise by the yellow ones. When the signal to noise ratio is higher than the predefined threshold the result is “good” and a green dot is shown below the corresponding frequency; otherwise, the result is “not good” and the shown dot is red. To obtain a “good” result, the signal level must exceed the threshold value defined in setup (-5, -10 or -15dB SPL), otherwise it is not considered and the signal indicator turns gray.

Beside the graph histograms show signal and noise.

Progress, Noise and SNR indicators are updated in real time while the test is in progress. At the end of the test the numerical result will be shown.

Press **SETUP** (see picture beside) to open the test setup page.

The menu allows for the setting of a few test parameters:

### **Adjustable parameters:**

#### **Time Out**

Maximum search time per single frequency. By touching this field the maximum waiting time, between 2 and 16 seconds, can be selected. The Auto mode is also available with a time varying from a minimum of 2 to a maximum of 8 seconds: the algorithm checks the signal-to-noise ratio in real time and when it is at least 1 dB higher than the set threshold value, it passes directly to the frequency subsequent speeding up the execution of the test.

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This field allows to select the pure tone signal intensity for F1.

Possible options:

L1 = L2, the level of F1 will be the same as the level set for F2

L1 = L2 + 10 dB, the level of F1 will be the same as the level set for F2 plus 10 dB

Auto: the level of F1 will be determined based on the following formula:  $L1 = 0.4 L2 + 39$  dB. When selecting this option, the stimulus level is automatically adjusted to obtain the maximum cochlear response.

#### **Stimulus Level F2**

This field allows to select the pure tone signal intensity for F2, between 30 and 70 dB SPL.

#### **Artefact Level**

In this field the maximum acceptable noise level can be selected. When the noise exceeds the set value, the measurement is rejected (artefact). By touching this field a value between 20 and 60 dB SPL can be selected.

#### **Frequency retest**

When enabled, it tests again the specific frequency, by slightly change the value, in case of "refer".

#### **Number of frequencies per pass**

Minimum number of frequency bands that must exceed the threshold SNR value set to obtain the test pass.

#### **Frequencies and SNR**

The value of the signal to noise ratio for each frequency can be set between 3, 6 or 9 dB, or the frequency band can be disabled by setting "-".

A maximum of 3 sets of predefined values can be saved by keeping the numbered fields on the right side of the screen pressed for a few seconds and until they turn green. In addition to the 3 user-definable protocols, there are 2 presets: 55/65 and 70/70. These are used for screening tests, perform 5-frequency checks at fixed levels, and cannot be edited

## 4.10 Otoscopy

The R16M audio-tymp unit can provide a very efficient tool to analyze the external ear conditions, when coupled with the RE-VO otoscope for the image acquisition of the external ear-canal and of the tympanic membrane. Like all the other audiological tests, images can be stored in the local database and later on transferred into the MDS Software.

### Connection

The otoscope is equipped with a dedicated cable: all you have to do in order to connect it to the device is inserting the mini-USB connector in the otoscope and the USB connector in the R15C USB Host; it automatically switches on when an examination is started, and you do not need to press any buttons to start.

### Speculum insertion

The otoscope is supplied with a set of multiple size disposable specula making it suitable for both young and grown-up patients.

Place the speculum on the otoscope tip and press it until it is fully inserted. In case the image appearing on the screen is not centered, softly move the speculum tip around, until the correct position is achieved.

Specula are single-use and need to be replaced each time after using.

Only specula supplied by MRS can be used.

In case any other part of the device comes into contact with the patient, just disinfect.


The unit is made of biocompatible materials.

### Test performance


In order to start the test, switch the device on and press the icon illustrated on the side, on the screen.

When the test is started all LEDs light up. Brightness is automatically adjusted to obtain an optimal image quality however; you can also manually change it at your will by using the + and – buttons positioned on the handle. By keeping the – button pressed for a few seconds you can switch from 2 LEDs mode to 4 LEDs mode and vice versa

Hold the otoscope with its focusing ring facing downwards so that the image appears straight (i.e. not rotated). The focus ring allows for both focus adjustments and image blowing-up.

Softly insert the speculum into the ear canal and adjust the focus using the focus ring until you obtain a sharp image.

For capturing and saving the image into the gallery, press the screen on the image or icon illustrated on the side. You can save up to 9 images for both the right ear and the left ear. In case you capture more images, the last 9 ones will be saved overwriting the previous ones.


To switch side, press the R or L icons.

To review the captured images, press the Gallery icon.

The gallery shows all saved images and in order to see them in extended view you just need to touch them. Once a full-screen image is obtained you can further blow it up by double-clicking on the screen. The zoomed image will become larger than the screen, but you can still see it all by keeping your finger firmly on the screen while rotating it in the different directions. To go back to the full-screen image, double-click on the screen again.

By pressing the **Gallery icon** you go back to all saved images in small size.

Press the **Otoscope icon** to leave the Gallery and go back to the otoscope main screen.

In order to delete an image keep your finger pressed on the small image you wish to delete for a few seconds.

### Test saving and exit

Pressing the **EXIT** button, a dialog box will pop up and ask if you wish to save and/or print the test out on the thermal printer.

Press the icon corresponding to your wish: if you confirm you want to save the test, the R16M will automatically address you to the patient database where the patient being tested can be selected and confirmed. For additional information on the patient database, please refer to section #6 in this manual.

## 5. Setup settings

After switching on the device, set the parameters that affect the test operation.


Press the **SETUP** icon (see adjacent image) to enter the setup page, there the settings are grouped as follow: GENERAL; IMPEDANCE; OAE.

### 5.1 General Setup

Press one of the tabs to access the respective Setup page:

## Bluetooth

Press the **BLUETOOTH** icon (see adjacent image) to activate database data transmission to PC, if set to OFF transmission is disabled.

## Autosave

Press the **AUTOSAVE** icon (see adjacent image) to enable the autosave mode. If enabled, once selected a patient, every performed test will be automatically stored into the database after pressing the **EXIT** icon; if disabled a pop up window will be displayed asking for confirmation.

## Language

Through the **LANGUAGE** button (see adjacent image) you can select the interface language. By pressing this icon a menu with all the available languages will be opened.

## Rotation block

Some of the tests provide for their visualization either horizontally or vertically: the visualization is automatically set by rotating the instrument. This option allows to block the rotation keeping the display horizontal.

## Wi-Fi

By pressing the **WIFI** icon (see adjacent image) a window allowing the connection to a Wi-Fi network will. This icon is active only when a Wi-Fi USB dongle is inserted in the USB host port. This option can be used to check the connection with the wireless printer.

## Power management

In order to improve the battery operating autonomy, you can set 3 different power modes: energy saving, balanced mode and performing mode. Based on the selected mode the display backlight and the unit automatic switching off will change. When the unit is connected to the power supply and the performing mode is selected, the display remains constantly backlit.

## PIN code change

By pressing this icon, you can change the PIN code by typing in any 4-digit numeric code or disable the PIN code request at power on. Please note that in order to change or disable the PIN you need to type in the one currently in use first.

## Patient data mandatory fields

Pressing this key gives access to a menu where it is possible to define the minimum fields necessary for the creation of a new patient. It is possible to set only a numeric code or activate the fields like name, surname, etc.

## Print-out


By pressing the **PRINT** icon (see adjacent image) you can set the print options for the optional portable printer.

You can select the printing-out of the following fields:

- Unit serial number and logo
- Patient's data
- Calibration expiry date
- Lines for comments after the heading
- Lines for comments after every test


## Report header customization

You can customize the report header you print by importing a graphic file that contains your logo or information about your company. The graphic file must be in BMP format, black and white (1 bit x pixel) and have a maximum size of 600x400 pixels. Create a folder named PrinterLogo in the root of the pen drive, copy the file or files to be used inside and insert it in the USB HOST port of the device. It is necessary to copy the file to the device memory to be able to use it for printing.

Pressing the **LOGO** button (see adjacent image) will open the page for selecting the file

Using the **ARROW** keys, it is possible to choose the graphic file from those present both in the device and in the pen drive

If the selected file is on the pen drive, the + key will turn orange, pressing it will copy it to the device memory.

Once the desired file has been selected, press the CONFIRM key (see adjacent image) to set it as the print header.

The DELETE button will allow you to delete the selected file.


If the printer is connected, using this key it is possible to print the test page with its logo.


## Licenses

Press the **LICENSES** icon (see adjacent image) to check the active test licenses. The following is displayed:


To add a new license, press the **ADD NEW** button and enter the license code:


The **ADD NEW** button should be pressed again to confirm.

The device recognizes the validity of the license and will display a confirmation message stating the test that has been added.

An error is displayed if the license is not valid.

If there are no active licenses, the tests can be used in demonstration mode up to 10 times.

The number in brackets shows the number of demonstrations remaining.

## Date and time

26  
Monday  
September

Press the **DATE** icon (see adjacent image) to set date and time.

Use the arrows to adjust the date and time.

It is possible to change the date format and related separator characters using the appropriate field using the left-right arrows

Press this icon to confirm.

## About


Pressing the **ABOUT** icon displays the devices' serial number, firmware and software versions. In addition, when entering this page, the device executes an internal test of the main hardware functionalities and shows the results.

## 5.2 Impedance Setup

Press the **IMPEDANCE** tab to enter the respective setup page:

### Sensitivity

Press the **SENSITIVITY** icon (see adjacent image) to change the sensitivity during reflexes test. When **HIGH** sensitivity is selected, reflexes with at least 30% of the points over the 50 ul threshold are detected. When **LOW** sensitivity is selected, reflexes with at least 30% of the points over the 25 ul threshold are detected.

### Autorun

Press the **AUTORUN** icon (see adjacent image) to activate the AUTORUN mode (Autorun ON). This feature enables the automatic test start. For correct use, follow these steps

- 1) Activate the Autorun function in the Impedance Setup menu
- 2) Enter the desired test
- 3) Wait until the device detects the Probe Open status
- 4) Insert the probe into the ear
- 5) Wait until the test starts

At the end of the test, in case you want to run another test, press the start button, or remove the probe. Let the software detect the Probe Open status and reinsert it into the ear to automatically start the new exam.

### Cavity Cal.

Press the **CAVITY CAL.** icon (see adjacent image) to enter the Compliance Calibration page. For further information refer to section 5.4 of this manual.

### Reflex slope

Through this icon (see adjacent image) you can activate the reflex rising and falling time while zooming it. The value is expressed in ml/mec

### Reflex display

By pressing the **VIEW** icon (see adjacent image) you can choose whether to visualize the test results either in positive or negative values.

## Y, B, G curves visualization

Press the YBG icon (see adjacent image) to change view in tympanometry 1000 Hz Admittance (Y); Conductance (G); Susceptance (B).

## Contralateral transducer

Press the **INSERT** icon (see adjacent image) to select the Contralateral transducer in use.

This function allows switching between different contralateral transducers. Ensure that any transducer to be connected is calibrated to the device. To change the settings keep the button pressed for some seconds.

## Tympanogram classification

When the tympanogram classification is enabled the device analyses the curve and provides for its automatic identification as A, B, C or AD based on Jerger classification.

## Calibration

Press the **CALIBRATION** icon (see adjacent image) to enter the Impedance calibration page.

**Only authorized personnel are allowed to perform calibration.**

## 5.3 Compliance Calibration

It is advisable to run the compliance calibration check daily, using the cavities supplied.


To perform the calibration check, enter either tympanometry or reflex assessment mode.

Using a yellow clinical ear tip, insert the probe into the test cavities and ensure that the compliance values read on the device correspond to the size of the cavity selected.


If measured values differ by more than 0.2 ml, it is recommended to follow this procedure to reset the compliance calibration:


- 1) Press the SETUP icon in the main menu and select impedance.
- 2) Select the CAVITY CAL. icon and enter the calibration menu. Follow the instructions displayed on the device.

- 3) Insert the probe into the cavity highlighted (5 ml)


- 4) Wait until the end of the control process (yellow indication on) without moving the probe.

- 5) The indicator will turn green at the end of the calibration, proceed to the next cavity as indicated on the screen.


- 6) After the calibration is completed, press the EXIT icon.

## 5.4 OAE Setup

Press the **OAE** tab to enter the respective setup page:

### End of test criterion after obtaining the PASS

With this key (see adjacent image) it is possible to establish the criterion for ending the test in case of Pass: it is possible to end the test just when the Pass is obtained or, even if the pass was obtained, continue to the end of time and see overall result.

### DP-OAE test frequency scan

Pressing this button (see adjacent image) changes the scanning order of the DP test frequencies: from the minor to the major or vice versa.

### OAE Probe check

This function allows you to check the operation of the OAE probe. Press the key to access the verification menu, it is necessary to insert the probe in the test cavity provided to start the verification.

The functionality will be checked and the result will be provided on the screen; if OK, the tests can be carried out, otherwise the probe must be checked for: tip cleaning; tip correctly inserted into probe body; front part probe coupling (both brackets hooked), see chap. 11.1 for more information.

If the malfunction persists, contact technical support.

## Extension of the maximum ear volume limit

The device automatically adjusts the levels of stimuli adapting to the volume of the ear, the maximum limit is about 2.8ml. In some cases, it may be necessary to perform the test even if a higher volume is detected (e.g. perforated eardrum). With this key it is possible to extend the test up to volumes of about 4ml; in the latter case the levels will be generated as if you were in the presence of a 2ml volume.

## Minimum OAE level in DP test

In the DP test the OAE must be higher than a minimum threshold level to be considered as valid and then proceed to evaluate the signal to noise ratio; this level can be set using this button at -5, -10 or -15dB SPL.


5.6 Audiometry, Impedance and OAE Calibration

## 5.5 Impedance and OAE Calibration

Entering the relevant setup menu (audiometry, impedance or OAE) and pressing the **CALIBRATION** icon (see adjacent image), you will be able to access to the calibration page.

**Only authorized personnel are allowed to perform a new calibration.**


The device calibration must be performed by trained technical staff only. Information concerning how to perform a calibration are available on the Technical Manual

## 6. Patient database management

The device has an internal memory which can store patient data and test results. The database is accessed by pressing the **PATIENT** button on the home screen.

In the database, all data stored can be displayed.

### Database

The left column lists patients, the center column shows test sessions associated to the patient, and the right-hand column displays demographics.

Using the **ARROW** keys to navigate between patients, results for the highlighted patient will be displayed

Patients list can be sorted by name or company name (if inserted in the patient data). Press the icon shown on the left to change the sorting option.

### Patient data

The right side of the screen displays the selected patient's personal data. The lower part contains some icons that activate the specific functions described below.

#### Create new patient.

Use the + button to add a new patient to the database. When executing this command, a virtual alphanumeric keyboard will appear which allows data to be entered.

To display the accented characters, simply press and hold the key for a few seconds: a pop-up will appear with the alternatives (e.g. hold down the A key to enter Å, Ä, Å, Ä, ...) Using the options key, you can select different keyboard layouts with the corresponding character sets (e.g. Qwerty, Arabic, Cyrillic ...)

By shortly pressing the **DELETE** icon you will delete the highlighted test. Once all tests have been deleted, the patient can be removed from the database by pressing the **DELETE** icon again. As an alternative you can delete a patient and all the related tests by keeping the **DELETE** key for a couple of seconds and then confirming the action.

Multiple test elimination: if you wish to delete multiple tests belonging to the same patient, select the patient, keep the first test you wish to delete pressed for a few seconds and until it turns orange, then add the other tests you wish to delete by touching them shortly. When all the selected tests will turn orange, you will be able to delete them all by pressing **DELETE**. The action can be aborted by pressing the Exit button.

**Patient data editing.** Pressing the **EDIT** icon allows data editing for the selected patient.

**Test results display.** Press the **VIEW** icon (see adjacent image) to display a specific test result for the selected patient.

**Patient activation.** Press the icon shown on the left to activate the selected patient. The selected patient's name will be displayed on the top left of the display once you enter one of the pages that allow the tests execution. Tests that have been performed will be saved to the associated patient.

## 7. Print management

The test results can be printed through the PC Suite Resonance MDS, or directly from the unit to the optional portable wireless printer.

To print out via the MDS, Transfer the data from the device to the PC Suite Resonance MDS and then follow the instructions related to the "Session report" MDS function.

If you wish to print directly from the unit, switch the printer on and insert the Wi-Fi dongle in the R16M USB port. Please note that the two devices require 20-30 seconds to connect.

Once unit and printer are connected, enter the patient database menu, select the test you wish to print out and press the **PRINT** icon you will see on the display. It is possible to print report during the text exit by pressing the print button on the pup up menu.

The wireless printer uses the 2.4 GHz Wi-Fi technology. Using it in areas where too many Wi-Fi networks, or other sources of electromagnetic interference such as mobile phones or other similar devices can be found, can provoke malfunctions.

## 8. Data transmission to PC

Data transfer from the device to PC is possible via Bluetooth connection.

To transfer test results stored within the R17A and subsequently view, edit, print or store them in a NOAH database, the PC suite Resonance MDS must be installed on to a computer. The installation file named **MRS Resonance Setup.exe** is supplied on the USB pen-drive which includes this operating manual.

Refer to the installation & operating instructions provided with the PC suite Resonance MDS for further details.

### How to associate R16M to PC via Bluetooth

1. Activate the Bluetooth connection on the device by pressing Bluetooth button present in the General Setup page (see section 5.1 of this manual).
2. Make sure the Bluetooth data transmission is active on your PC.
3. Open **Windows settings** on PC.
4. Select the **Devices** icon.

5. Open **Bluetooth** devices management feature.
6. Wait until your PC detects the device. This may take several minutes. When finished, an icon showing the serial number of the device to be associated will be displayed.
7. Check that the serial number displayed corresponds to the one on the device identification label. If it matches, press the icon on screen to complete the association.
8. The unit will generate a **PASSCODE** which you will see displayed on the screen.
9. The **PASSCODE** will be displayed on the PC screen, too.
10. If the two codes are the same, press the confirm button on the unit first and then on the PC.
11. When this process is finalized, the unit is associated.

As an alternative, data can be transferred via cable through the unit mini USB port. If you choose this option be aware that by connecting a Personal Computer to a Medical Device you create a system which, if used in a patients' area need to comply with all safety requirements. In case the R16M is connected to other devices equipped with their own external power supplies, the compliance of the whole system to the IEC 60601-1 standard needs to be verified under the technician's connecting the units responsibility. In case insulation devices are used, they have to comply with the IEC 60601-1 standard. Each person connecting an R16M to other external units creates a system and is held responsible for the compliance of the system itself to the IEC 60601-1 standard.

# 9. Main Labels

## Device labels

TOP HOUSING

RIGHT SIDE

FRONT VIEW


LEFT SIDE

BOTTOM HOUSING

## 10. Messages and warnings

Save the test in the database?  
Is it possible to perform the print out. The button is active only if the printer is active and connected to the Wi-Fi adapter.

Save current settings in setup?


Update also the calibration date?


Calibration expired. Recalibration required. Continue anyway?

Instruct the patient! For the ETF test it is necessary to instruct the patient to swallow.


Test ended incorrectly. It was not possible to reach the set pressure for tubal function testing. Check the tip seal.

Attention, probe open. Check the probe

Attention, probe closed. Check the probe


Delete the exam from database?

Delete the patient from database?


Attention, optional test license is not correct


Optional test license is correct.


Attention, battery is running low. Connect the device to its medical power supply to recharge and continue the test.


Delete the selected point?


The noise level is too high for the test to be performed. Check the probe correct positioning and/or the selected noise level in the Setup menu.


The connected probe has no serial number. You need to confirm whether you want it to be permanently coupled to this specific device.


An error occurred while reading the probe data. Check the probe connector.

The connected probe is not suitable for the selected test. The probe needs to be replaced with the correct one.


No probe is connected. Connect a probe.


The probe is not properly inserted into the ear. Check the tip size.


The connected probe is coupled to another device. The test can be run anyway


The internal probe memory cannot be updated with the calibration data! Exit without saving?


Pair the connected probe with this device? This operation will configure the probe to work permanently with the device

Please check the probe! Check that transparent probe tip is clean from any impurity and both plastic brackets are hooked. Refer to section 11.2

Probe does not work!

The free memory for saving exams is running out. It is possible to free up memory by transferring the data to the MDS suite and deleting them from the device

# 11. Care and maintenance

Switch off the device and disconnect the power plug before cleaning!

Always disconnect the mains conductor during the cleaning process, and be careful that no fluid enters inside of the device or the accessories; no alcohol or spirits should be used.

The device requires cleaning the housing with a damp cloth (do not use aggressive cleaning agents); avoid liquid ingress.

During testing operations, headset pads and probe tips come into direct contact with the patient. Although said components are biocompatible, they are not sterile and they must therefore be subject to special care:

After each test, it is necessary to disinfect the headset cushion before use on a new patient. This includes physically cleaning using of a recognized disinfectant. Individual manufacturer's instruction should be followed for use of this disinfecting agent to provide an appropriated level of cleanliness.

The insertion headsets tips are disposable: use only once and substitute them before use with a new patient.

Calibration and any possible repair, installation and/or update of the device must be strictly carried out only by specialized technical personnel, duly authorized by the R15C manufacturer and in full compliance with the terms and conditions specified in the Technical Manual.

Standard maintenance and device calibration operations are to be regularly performed at one-year intervals. Standard maintenance operations cover the inspection of the transducer cable conditions, of the mains power supply cable, of the earphone pads as well as an inspection check to ensure the device's outer section and/or cabinet are in good working order. All device calibration operations must be strictly performed only by qualified technical personnel, duly authorized by the manufacturer and implementing the appropriate apparatus for checking the levels and frequencies of all the output signals by the transducers supplied with the device. Once device calibration testing using the special purpose apparatus is complete, it is necessary to run a device function test.

Should transducers require replacement, use only parts provided by Resonance. Upon replacement of transducers, the entire system will require re-calibration

The device is provided with a lithium-ion battery type 1SP NCR-18650A. Said battery is fitted into the device and possible replacement thereof must strictly be performed only by qualified technical personnel. In case of replacement, use only the same battery type.


## 11.1 Impedance Probe maintenance

Ensure the probe is clean prior to use or calibration. Cerumen may have entered the probe tip and should be removed. To remove wax, the conical end to the probe may be unscrewed. Cerumen may also have entered the probe tubes. To clean the tubes, carefully remove them from the probe body. Do not pull the tubes from the adaptor, hold the adaptor base firmly (as shown in the diagram).

In this case, unscrew the conical part of the probe and remove any earwax from it.


Earwax may also have penetrated in the probe tubes which are visible after removing the conical part. The two tubes administer the probe signals.


In order to clean the tubes, remove them from the probe body by pulling them out (please note that the two tubes are fixed on a base adapter which will be removed with them). Do not remove the adapter pulling the tubes. Use your fingers to pull the adapter.


**Wrong way** to extract the inside part of the probe

**Correct way** to extract the inside part of the probe

Use the pipe cleaner provided to clean the two tubes by simply sliding it into them.



To reassemble the probe, put the tubes with the adapter into the probe body and then screw the conical part onto the probe body.

Place the adapter into the probe body by matching the two parts.

Care should be taken when reassembling the conical probe tip, it can easily be cross-threaded. Tighten to achieve an airtight seal.

When closed, verify that the two internal tubes are a little under the edge of the probe output as shown in the image below.

To ensure that the device is functioning correctly after reassembling the probe, a compliance calibration should be run.


## 11.2 Otoacoustic emission probe maintenance

Before using the probe, make sure it is clean and its transparent plastic tip is totally free from ear wax residues to prevent blockages.

To clean the probe proceed as follows:

1. Remove the front part of the probe by slightly widening the two plastic clips placed on the sides.

2. Remove the front part of the probe by slightly widening the two plastic clips placed on the sides.


3. Pull out the transparent part of the tip.


4. Clean the 3 tip channels using the supplied cleaning tool.

5. Make sure the 3 channels are completely free from any kind of debris.

6. Place the transparent probe tip back making sure the recess on the tip is properly matched with the corresponding tooth on the front part of the probe. Please note that in case the tip is inaccurately positioned, closing the probe is not possible.

7. Place the front part of the probe tip back into place making sure it is perfectly coupled to its body and both plastic brackets are hooked.

## 12. Safety regulations

### 12.1 Electrical Safety

The device applied parts comply to type BF of IEC 60601-1, in any case, the device is not intended to be used in environments dealing with explosive material or equipment.

### 12.2 Device handling and operation

The device should be checked once daily, as per operations described on section 11.

The device is designed to be handled and operated by trained personnel such as Audiologist, ENT doctors or personnel with similar qualifications ONLY.

### 12.3 Measurement Safety

As a manufacturer, we recommend to perform an annual inspection and calibration in order to ensure safety and quality of the measurement.

Annual calibration and maintenance should only be carried out by RESONANCE authorized personnel. M.R.S. will not be liable for any failure to comply with the specified inspection date, according to the Medical product law. The use of uncalibrated and uninspected device is strictly prohibited

## 13. Technical specifications

| Model                                | R16M       |            | R16M-D     | R16M OAE                                     | R16M-D OAE                                   |
|--------------------------------------|------------|------------|------------|--|--|
| Code                                 | 4300104100 | 4300104145 | 4300104110 | 4300104115<br>4300104125<br>4300104135       | 4300104120<br>4300104130<br>4300104140       |
| <b>Test</b>                          |            |            |            |  |  |
| Tympanometry (226 Hz)                | ●          | ●          | ●          | ●  | ●  |
| Tympanometry HF (1000 Hz)            | ●          | ●          | ●          | ●  | ●  |
| Reflexes Ipsi                        | ●          | ●          | ●          | ●  | ●  |
| Reflexes contra                      | ●          | -          | ●          | ●  | ●  |
| Quick                                | ●          | ●          | ●          | ○  | ●  |
| Decay                                | ○          | -          | ●          | ○  | ●  |
| Decay Free Run                       | ○          | -          | ●          | ○  | ●  |
| ETF intact                           | ○          | -          | ●          | ○  | ●  |
| TEOAE<br>Screening/Diagnostic        | -          | -          | -          | ● 4300104125<br>○ 4300104115<br>● 4300104135 | ● 4300104130<br>○ 4300104120<br>● 4300104140 |
| DPOAE                                | -          | -          | -          | ○ 4300104125<br>● 4300104115<br>● 4300104135 | ○ 4300104130<br>● 4300104120<br>● 4300104140 |
| <b>Accessories</b>                   |            |            |            |  |  |
| Impedance Probe*                     | ●          | ●          | ●          | ●  | ●  |
| Otoacoustic emission probe*          | -          | -          | -          | ●  | ●  |
| Insert for Contra*                   | ●          | -          | -          | ●  | -  |
| IP30 Insert earphones for<br>Contra* | ○          | -          | ●          | ○  | ●  |
| Impedance Kit tips*                  | ●          | ●          | ●          | ●  | ●  |
| OAE Kit tips* and test cavity        | -          | -          | -          | ●  | ●  |
| Calibration cavity                   | ●          | ●          | ●          | ●  | ●  |
| USB pen-drive with User<br>Manual    | ●          | ●          | ●          | ●  | ●  |
| Medical CE approved power<br>supply  | ●          | ●          | ●          | ●  | ●  |
| Cradle                               | ●          | ○          | ●          | ●  | ●  |
| MDS software                         | ●          | ●          | ●          | ●  | ●  |
| External USB keyboard                | ○          | ○          | ○          | ○  | ○  |
| Internal Bluetooth module<br>2.1/4.0 | ●          | ●          | ●          | ●  | ●  |
| USB Cable                            | ○          | ○          | ○          | ○  | ○  |
| Wireless Printer                     | ○          | ○          | ○          | ○  | ○  |
| RE-VO Otoscope                       | ○          | -          | ○          | ○  | ○  |

● = standard ○ = optional - = not available \* = applied parts

Applied parts are constituted by the transducers that convey the test signals. Said transducers are applied so that they go into direct contact with the patient. In the foregoing list, they are highlighted in bold script and are counter-marked with an asterisk (\*).

|                     |   |
|---------------------|---|
| Reference standard: | Impedance: IEC 60645-5; ANSI S3.39 Type 2<br>Otoacoustic emission: EN 60645-6 Type 2  |
| Safety:             | IEC 60601-1 (2006) + A1 (2012) Class II Type BF   |
| EMC:                | IEC 60601-1-2 (2014)  |
| Medical CE Marks:   | Class IIa<br>The Manufacturer's Quality Management System complies to the Annex II of Medical Device Directive 93/42/EEC. The notified body IMQ, Identification N. 0051, has approved the Quality System. |

| Impedance   |  |                       |                                  |                     |
|---|--|-----------------------|----------------------------------|---------------------|
| Tympanometry (automatic, absolute and compensated)                  |  |                       |                                  |                     |
| Multifrequency Tympanometry (Y-B-G track visualization)             |  |                       |                                  |                     |
| Quick test (performing tympanometry and reflexes in rapid sequence) |  |                       |                                  |                     |
| Diagnostic Test: Decav, ETF intact, Free run (depending on model)   |  |                       |                                  |                     |
| Tympanometry  | Probe tone frequency: 226 Hz, 1000 Hz, 1000 Hz accuracy $\pm 1\%$                            |                       |                                  |                     |
|   | Prob tone intensity: 85 dB SPL 85 dB SPL $\pm 2$ dB (226 Hz); 75 dB SPL $\pm 2$ dB (1000 Hz) |                       |                                  |                     |
|   | Pressure range:  |                       |                                  |                     |
|   | Pressure accuracy: $\pm 10$ daPa oppure $\pm 10\%$   |                       |                                  |                     |
|   | Pressure generation speed: 200 daPa/s, and automatic   |                       |                                  |                     |
|   | Compliance range: 0.02 ml – 7 ml   |                       |                                  |                     |
| Compliance accuracy: $\pm 0.1$ ml or 5%                             |  |                       |                                  |                     |
| Reflex and Decay test max frequencies and max levels                |  |                       |                                  |                     |
| Tone Freq. Hz   | IPSI (dB HL)   | Contra Insert (dB HL) | Contra Insert phone IP30 (dB HL) | Contra DD45 (dB HL) |
| 500   | 100  | 110                   | 110                              | 105                 |
| 1000  | 100  | 110                   | 110                              | 110                 |
| 2000  | 100  | 110                   | 110                              | 110                 |
| 4000  | 100  | 100                   | 110                              | 110                 |
| Reflexes and Decay test   | Stimulus duration  |                       |                                  |                     |
| Tone frequency accuracy: $\pm 1\%$                                  | Reflexes: 0.5 sec, 1 sec   |                       |                                  |                     |
|   | Decay and Decay free run: 10 sec, 20 sec   |                       |                                  |                     |

Technical descriptions more detailed and complete as wiring diagrams, parts list, instructions for calibration and other materials necessary for the repair and maintenance of the instrument, are contained in the service manual which will be provided, upon request, technical personnel authorized by the manufacturer to carry out such operations

| Otoacoustic Emission  |                              |
|-----------------------|------------------------------|
| Number of probe:      | 1                            |
| Acquisition           | 24 bit 44.1 Ksps             |
| TEOAE                 |                              |
| Stimulus level:       | 70 ÷ 90 dB peSPL             |
| Stimulus frequency:   | 66 Hz                        |
| Stimulus type:        | Non- linear Click            |
| Measurement interval: | 512 samples, da 5 a 15 msec. |
| Bandwidth:            | from 750 Hz to 5KHz          |
| Artefact noise level: | from 20 to 60 dB SPL         |
| Protocols:            | 3, user's definable          |

| TEOAE Screening          |  |
|--------------------------|--|
| Pass stop criterion:     | Statistical analysis of the signal polarity                          |
| Max test time            | from 500 to 1250 stimuli (from 30 to 75 sec)                         |
| TEOAE Diagnostic         |  |
| Pass stop criterion:     | signal/noise ratio 3, 6 or 9 dB definable for 3, 4 or 5 frequencies. |
| Frequenze di indagine    | 1, 1.5, 2, 3, 4 Khz  |
| DPOAE                    |  |
| F2/F1 frequency ratio    | 1.22   |
| Reference frequency:     | F2, from 1.5 to 12 kHz, 15 frequencies Max, user's definable         |
| Stimulus level L2:       | from 35 to 70 dB SPL step size 5 dB                                  |
| Stimulus level L1:       | User's definible: automatic (L1=0,4* L2+39 dB), L1=L2, L1=L2+ 10 dB  |
| Artefact noise level:    | from 20 to 60 dB SPL   |
| Detection criterion:     | Signal noise ratio 3, 6 or 9 dB frequency definable                  |
| Stop criterion:          | pass frequency number: from 2 to the number of enabled frequencies   |
| Test time per frequency: | automatic or user definable, 2 to 16 sec                             |
| Protocols:               | 2 fixed (70 dB and 65-55 dB), 3 user's definable                     |

Technical descriptions more detailed and complete as wiring diagrams, parts list, instructions for calibration and other materials necessary for the repair and maintenance of the instrument, are contained in the service manual which will be provided, upon request, technical personnel authorized by the manufacturer to carry out such operations

**Dimensions**

W x D x H: 180 x 135 x 41 mm

**Weight:** 600 g net**Environmental conditions****Transport and storage:**

Temperature: from -20 to \*50°C  
Humidity: from 20 to 90% without condensation  
Pressure: from 500 hPa to 1060 hPa  
Temperature: from 20 to 90% without condensation  
Pressure: from 700 hPa to 1060 hPa

**Power Batteries:**

Internal rechargeable Lithium-Ion battery 3,6V-3,1 Ah  
Battery status indication

**Battery duration:**

Normal use, up to 8 h  
Continuous use (test performance): 3 h

**Battery charge time:**

Less than 4 h; 80% in 2 h

**External power supply:**

Medical power supply input 100/240 VAC 50-60 Hz- output +5 VDC/2A

**Construction**

Self-extinguishing plastic ABS cabinet

## 14. Electromagnetic compatibility

| Guidance and manufacturer's declaration– Electromagnetic emissions  |            |  |
|---|------------|--|
| L' R16M is suitable for use in the specified electromagnetic environment. The purchaser or user of the R16M should assure that it is used in an electromagnetic environment as described below: |            |  |
| Emissions test  | Compliance | Electromagnetic environment  |
| RF emissions<br>CISPR11   | Group 1    | This R16M uses RF energy only for its internal function. Therefore, the RF emission is very low and not likely to cause any interference in nearby electronic equipment.   |
| RF emissions<br>CISPR11   | Class B    | This R16M is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions<br>IEC/EN 61000-3-2  | Complies   |  |
| Voltage fluctuations/flicker emissions<br>IEC 61000-3-3   | Complies   |  |

| Guidance and manufacturer's declaration– Electromagnetic immunity   |   |                             |  |
|---|---|-----------------------------|--|
| This R16M is suitable for use in the specified electromagnetic environment. The purchaser or user of the R16M should assure that it is used in an electromagnetic environment as described below: |   |                             |  |
| Immunity test   | IEC 60601-1-2<br>Test level   | Compliance<br>level         | Electromagnetic Environment  |
| Electrostatic<br>discharge (ESD)<br>IEC/EN 61000-4-2  | 8 kV contact<br>2/4/8/15 kV air                                     | IEC 60601-1-2<br>Test level | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30% |
| Radiated<br>electromagnetic<br>field<br>IEC 61000-4-3   | 3 V7m<br>80 MHz to 2.7 GHz  | IEC 60601-1-2<br>Test level | Portable and mobile FR communications equipment should be used no closer to any part of the R15C including cables. Minimum distance 30 cm.   |
| Electrical fast<br>transient/burst<br>IEC/EN 61000-4-4  | 2kV for power supply<br>lines<br>1 kV for input/output<br>lines >3m | IEC 60601-1-2<br>Test level | Mains power quality should be that of a typical commercial or hospital environment.  |
| Surge IEC 61000-4-5   | 0.5/1 kV differential<br>mode<br>0.5/1/2 kV common<br>mode          | IEC 60601-1-2<br>Test level | Mains power quality should be that of a typical commercial or hospital environment.  |
| Conducted<br>disturbances<br>induced by RF fields<br>IEC 61000-4-6  | 3V<br>150 kHz to 80 MHz   | IEC 60601-1-2<br>Test level | Portable and mobile RF communications equipment should be used no closer to any part of the R15C including cables. Minimum distance 30 cm.   |

|   |  |                                      |  |
|---|--|--------------------------------------|--|
| <p>Voltage dips short interruptions and voltage variations on power supply input lines IEC 61000-4-11</p> | <p>0% Un. For 0.5 cycle<br/> 0% Un. For 1 cycle<br/> 70% Un. For 25 cycles<br/> 0% Un. For 5 s</p> | <p>IEC 60601-1-2<br/> Test level</p> | <p>Mains power quality should be that for a typical commercial or hospital environment. If the user of the R16M requires continued operation during power mains interruptions, it is recommended that the R16M be powered from an uninterruptible power supply or a battery.</p> |
| <p>Power frequency (50/60 Hz)<br/> Magnetic field<br/> IEC 61000-4-8</p>                                  | <p>30 A/m</p>  | <p>IEC 60601-1-2<br/> Test level</p> | <p>Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.</p>   |

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