


# DM-50 Cartridge Impression Material Comparison Chart

DM -50 CARTRIDGE 								
Material	Item Number	Cartridge Type	Color	Manufacturer	Package Contents	Viscosity	Setting Time	Shore-A Value
Flextime	JKR96B	DM-50	Pink / Violet	JKR Labs	8 cartridges, no tips	Low	2:30 min	35
SiliClone	10104	DM-50	Pink	Westone	8 cartridges + 36 MTL mix tips	Low	4-7 min	25
SiliClone	408	DM-50	Pink	Westone	8 cartridges, no tips	Low	4-7 min	25
Audiologist's Choice	01602	DM-50	Teal	Oaktree Products	8 cartridges + 36 MT4 mix tips	Medium	4m +/- 30s	25 +/- 2
Audiologist's Choice	01602NT	DM-50	Teal	Oaktree Products	8 cartridges, no tips	Medium	4m +/- 30s	25 +/- 2
Blue Ultra	BULTRA	DM-50	Blue	JKR Labs	8 cartridges, no tips	Medium	3m +/- 30s	37
Matrics	MATRR	DM-50	Teal	Matrics	8 cartridges + 36 MT4 mix tips	Medium	4m +/- 30s	25 +/- 2
Matrics	MATRRNT	DM-50	Teal	Matrics	8 cartridges, no tips	Medium	4m +/- 30s	25 +/- 2
Purple Stuff	PSR	DM-50	Blue	Matrics	8 cartridges + 36 MT4 mix tips	Medium	4-5 min	28
Purple Stuff	PSRNT	DM-50	Blue	Matrics	8 cartridges, no tips	Medium	4-5 min	28
Yellow Stuff IV	YS	DM-50	Yellow	Mid States Laboratories	8 cartridges + 32 MTYS mix tips	N/A	3-5 min	30
Yellow Stuff IV	YSNT	DM-50	Yellow	Mid States Laboratories	8 cartridges, no tips	N/A	3-5 min	30