


# DS-50 Cartridge Impression Material Comparison Chart

DS-50 CARTRIDGE 								
Material	Item Number	Cartridge Type	Color	Manufacturer	Package Contents	Viscosity	Setting Time	Shore-A Value
Silhouette Plus	7077790	DS-50	Teal	Siemens	8 cartridges + 36 S-50 MTM mix tips	Low	4-5 min	26 +/- 1
Silhouette Plus	7077741	DS-50	Teal	Siemens	8 cartridges, no tips	Low	4-5 min	26 +/- 1
SiliClone	10154	DS-50	Pink	Westone	8 cartridges + 36 S-50 MTM mix tips	Low	4-7 min	25
SiliClone	10153	DS-50	Pink	Westone	8 cartridges, no tips	Low	4-7 min	25
Accu-Flex II	3295	DS-50	Green	Hal Hen	8 cartridges, no tips	Medium	4m +/- 30s	35
Audiologist's Choice	01603S50	DS-50	Blue	Oaktree Products	8 cartridges + 36 S50 MTS mix tips	Medium	4m +/- 30s	25 +/- 2
Audiologist's Choice	01603S50NT	DS-50	Blue	Oaktree Products	8 cartridges, no tips	Medium	4m +/- 30s	25 +/- 2
Matrics NXT	M0836NXT	DS-50	Teal	Matrics	8 cartridges + 36 S-50 MTS mix tips	Medium	4m +/- 30s	25 +/- 2
Matrics NXT	M0836NXTNT	DS-50	Teal	Matrics	8 cartridges, no tips	Medium	4m +/- 30s	25 +/- 2
SiliClone Firm	10184	DS-50	Pink	Westone	8 cartridges + 36 S-50 MTL mix tips	Medium	4-7 min	25-30
SiliClone Firm	10183	DS-50	Pink	Westone	8 cartridges, no tips	Medium	4-7 min	25-30
Silicone X-ACT	10114	DS-50	Pink	Westone	8 cartridges + 36 S-50 MTM mix tips	Medium	3-4 min	20-22
Silicone X-ACT	10113	DS-50	Pink	Westone	8 cartridges, no tips	Medium	3-4 min	20-22
Matrics MB-34NXT	MB-34NXT	DS-50	Orange	Matrics	8 cartridges + 36 S-50 MTS mix tips	Medium	4m +/- 30s	34
Matrics MB-34NXT	MB-34NXTNT	DS-50	Orange	Matrics	8 cartridges, no tips	Medium	4m +/- 30s	34