

iQDry™

HEARING AID DRYER

Specifications



COMPONENTS

1. LCD display
2. Power button
3. USB Type-C port
4. Air inlet / outlet vent
5. Silicone feet



ICON INFORMATION

60-100%	40-59%	39-0%
Indicates humidity inside is DAMP	Indicates humidity inside is COMFORTABLE	Indicates humidity inside is DRY

DESCRIPTION

- ✦ **Shortest drying cycle of any dryer on the market. Only takes 15 minutes!**
- ✦ **Humidity indicator** shows the progress and result in removing moisture.
- ✦ **Ideal for rechargeable hearing aids.** Regulated temperature and extremely short drying cycle mean you can dry your rechargeable hearing aids in iQDry. Temperature max.: 104°F
- ✦ **Countdown timer with automatic shut off** tells you when the cycle is complete.
- ✦ **Stylish, compact design** easily fits on nightstand or countertop.
- ✦ **Built-in rechargeable battery so one charge can carry you through seven drying cycles.** Battery capacity reminder; USB Type-C charging compatibility
- ✦ **Easily accommodates two Behind-The-Ear (BTE) style hearing aids.**



Drying cycle time	15 minutes
Temperature	104°F maximum
Input / Output	5 V = 1 A / 5 V = 1 A
Rated capacity	800 mAh 2.96 W
Weight	3.8 oz.
Dimensions	2.8" x 3.3" x 1.8"



EAR TECHNOLOGY CORPORATION®

1-800-327-8547 | info@eartech.com | eartech.com