

**HEARING
INSTRUMENTS
ARE EVOLVING.
WE BELIEVE
INSTRUMENT CARE
SHOULD, TOO.**



I couldn't wear my hearing instruments very long because my hair sweats so bad. One didn't work at all and the other one wasn't very clear. After Redux did the drying, the one that didn't work started working and everything was a lot sharper, a lot clearer. Yeah, it made a heck of a difference.



CONTACT REDUX

Talk to your hearing care provider about how the proven effectiveness of Redux can save you money and improve the function of your hearing instruments.

REDUX CARE PROGRAM MEMBER

**EMAIL US AT HEARING@REDUX.COM
WWW.REDUX.COM**



FOR HEARING INSTRUMENTS



**HEARING
INSTRUMENT CARE
REINVENTED.**

www.redux.com | hearing@redux.com

OUR CARE PROGRAM

WHO BENEFITS?

Anyone using a hearing instrument. All devices have a level of moisture and bacteria that can impact functionality and lifespan. Physical activity or underlying health conditions such as obesity, diabetes, high blood pressure, or high cholesterol increase risk of water damage.

WHAT IS INCLUDED?

Choosing the Redux Care Program provides patients access to an innovative tool that can improve the functionality and lifespan of hearing instruments.

WHEN IS IT USED?


Anytime you have a few spare minutes to stop by your local hearing care provider.

WHERE IS IT FOUND?

Redux partners with a growing network of hearing care professionals to bring this Care Program to patients.




THE REDUX DIFFERENCE




FASTER

Redux is proven to dry hearing instruments in about 15 minutes with the push of a button.



DRYER

Redux applies endpoint drying technology to maximize moisture removal from treated devices.



SMARTER

Redux combines vacuum drying with humidity level measurement to vaporize liquid at a low temperature that is safe for the device.



VERIFIED

Redux provides real-time measure of moisture removed during the drying process and reports the total amount upon completion.