

# Tub & Accelerator Impression Material Comparison Chart

TUB & ACCELERATOR								
Material	Item Number	Color	Manufacturer	Package Contents	Number of Impressions	Viscosity	Setting Time	Shore-A Value
XL 100	XL100	Purple	JKR Labs	Tub & Accelerator	60	Medium	5-6 min	25
XL 200	XL200	Beige	JKR Labs	Tub & Accelerator	60	Medium	5-6 min	30
XL 300	XL300	White	JKR Labs	Tub & Accelerator	60	Medium	5-6 min	25
Blue Silicast	203A	Blue	Westone	Tub & Accelerator	50	High	5-8 min	30-32